


ITALIAN COMPANY

INDEX

06	POST-SURGICAL BRAS
18	BREAST AUGMENTATION BELT
22	LIPOSUCTION GARMENTS
32	BBL GARMENTS
38	TUMMY TUCK GARMENTS
48	REVÉE SLEEVES
54	FACE GARMENTS
58	REVÉE PADS
62	ATITHROMBOSIS STOCKINGS
66	MALE GARMENTS COLLECTION
72	REVÉE CARDIO
80	REVÉE SCAR GEL
84	REVÉE SUNCARE 50+
86	REVÉE PARTNER
88	REVÉE SUPPORT
94	REVÉE NEWS
96	DEFEDENCES











Company with **certified quality** management
system **TÜV Italia Srl No. 50 100 16173 Rev.001**



POST-SURGICAL BRA

Our Revée Bra Collection offers a full solution for post-operative recovery following breast augmentation, breast reduction, breast reconstruction.


It combines an exclusive design, engineered with the support of world reknown plastic surgeons with the highest-standards manufacturing process. Thanks to the unique characteristics of the Lycra® fabric (Latex free) it ensures optimal compression while providing maximum comfort to facillitate healing and repair of the tissues and infra-mammary structures (1). Patients who have undergone breast augmentation surgery have more success and satisfaction with their procedure when they wear a quality compression post surgical bra.













Ultra flat seam



No underwires



Breathable and quick to dry



3 front hooks



4 adjustable hooks

LAETITIA BRA

REV.0301



- 1 4 adjustable hooks on the shoulder straps.
- 2 Air Feel fabric on the back.
- 3 3 adjustable hooks on the front.
- 4 Extra breast support in Lycra[®].
- **5** Soft jaquard elastic bend by Revée[®] of 4.5 cm.
- **6** External seam to avoid irritating the nipple and areola.
- 7 Bra cups self-mold and perfectly adapt themselves in case of breast asymmetry.



Product code: REV.0301

Size: XS - S - M - L - XL - XXL - 3XL - 4XL

Colors:

Composition: 74% polyamide,

26% elastan

Washing: $\begin{tabular}{ll} \begin{tabular}{ll} \begin{tabular}$



SIZE TABLE AT PAGE 16






















Dry fast



Expandable cups



No underwires



Ultra flat seam



Front zip



4 adjustable hooks







ELODIE BRA

REV.0302



- 1 Closure: zipper and hooks.
- 2 4 positions adjustable hook straps.
- 3 Additional X-shape support in Lycra[®] on the back.
- 4 Soft jaquard elastic bend by Revée[®] of 4.5 cm.
- **5** Bra cups self-mold and perfectly adapt themselves in case of breast asymmetry.



Product code: REV.0302

Size: XS - S - M - L - XL - XXL - 3XL - 4XL

Colors: S9% polyamide,

11% elastan

Washing: ₩ ≥ ₩ ₩ ₩




SIZE TABLE AT PAGE 16















JULIE BRA

REV.0303



- 1 3 adjustable hooks on the front.
- 2 4 positions adjustable hook straps.
- 3 Additional X-shape support in Lycra[®] on the back.
- 4 Soft jaquard elastic bend by Revée[®] of 4.5 cm.
- **5** Bra cups self-mold and perfectly adapt themselves in case of breast asymmetry.



Product code: REV.0303

Size: XS - S - M - L - XL - XXL - 3XL - 4XL

Colors:

Composition: 89% polyamide,

11% elastan

Washing: Was



SIZE TABLE AT PAGE 16









LAETITIA

Size	Under breast
XS	66 - 70
S	71 - 75
М	76 - 80
L	81 - 85
XL	86 - 90
XXL	91 - 95
3XL	96 - 100
4XL	≥ 101

Under

breast
67 - 72
73 - 77
78 - 82
83 - 87
88 - 92
93 - 97
98 - 102
≥ 103

Under breast
67 - 72
73 - 77
78 - 82
83 - 87
88 - 92
93 - 97
98 - 102
≥ 103

*NOTES: MEASUREMENTS ARE INDICATIVE AND EXPRESSED IN CENTIMETERS

↑ Warnings:

- In order to guarantee its efficacy, tolerability and proper use, it should be worn with extreme care.
- In case of concern on the application procedure, contact the shop staff or the treating physician.
- specialist shop or threating physician recommends.
- The fabrics used are hypoallergenic, non-cytotoxic, and latex-free. is however advisable to check the recommended.

- components listed on the label to rule out any allergy to the fabrics used.
- In case a surgical wound is present, do not put into direct contact until complete healing.
- Close all the velcro and hook ■ Choose the correct size that your **fasteners** before washing in order not to damage other items and to protect the shaping feature of the product.

Tip: For a better post-surgical recovery, acquiring two bras of the right size is







REVÉE BREAST AUGMENTATION BELT



Revée Belt is a comfortable elastic band with an adjustable closure velcro to stabilize and avoid prosthesis mis-positioning. It is recommended after breast augmentation surgery in combinantion with Revee Bra.

Application fields:

Recommended after breast augmentation surgery or implants substitution



REV.0201



REV.0202









Enveloping belt



Designed by surgeons



Latex free



Adjustable velcro closure

REV.0201



- 1 Adjustable velcro closure
- 2 Latex free
- 3 8-cm-high hook-and-loop containment elastic by Revée®

Product code: REV.0201
Size: One-size-fits-all
Color:

Composition: 86% polyester,

14% rubber

Washing: 🕷 🛛 🖄

REV.0202



- 1 Adjustable velcro closure
- 2 Latex free
- 4 Contoured under the armpits

Product code: REV.0202

Size: S - M - L Color:

Composition: 87% polyamide,

13% elastan

Washing: 🖾 🔀 🖄 🛇

Size	Torax
S	70 - 85
M	86 - 105
L	> 105

*NOTES: MEASUREMENTS ARE INDICATIVE AND EXPRESSED IN CENTIMETERS










LIPOSUCTION GARMENTS

Lipo Revée[®] is a collection of multifunctional garments that meet all the specific requirements of an effective compression treatment, after liposuction interventions in one or more areas of the body: abdominoplasty, thighplasty, reshaping, torsoplasty, body lipofilling, and buttock lift.


Thanks to the Lycra® fabric, these garments are able to self-mold ensuring the best fit and promoting the healing process.















LIPO REVÉE

REV.0101



- 1 Triple front hook closure
- **2** Soft jaquard elastic bend by Revée[®] of 4.5 cm.
- 3 Abdominal support
- 4 Wide hygienic opening

Product code: REV.0101 - 0102 - 0103

Size: XS - S - M - L - XL - XXL

Color:

Composition: 74% polyamide,

26% elastan

Washing: 🚾 🔀 🕱 🛇











REV.0102



REV.0103









hooks





Best seller

XS 47 - 59 76 - 92 40 - 50 34 - 42 S 57 - 69 84 - 100 47 - 57 38 - 46 M 67 - 79 92 - 108 54 - 64 42 - 50 L 77 - 89 100 - 116 61 - 71 46 - 54 XL 108 - 124 68 - 78 50 - 58 87 - 99 XXL 97 - 109 116 - 132 75 - 85 54 - 62 3XL 107 - 119 124 - 140 82 - 92 58 - 66

*NOTES: MEASUREMENTS ARE INDICATIVE AND EXPRESSED IN CENTIMETERS



LIPO REVÉE

REV.0104



- 1 Triple front hook closure
- 2 Adjustable and removable straps
- 3 Abdominal support
- 4 Wide hygienic opening

Product code: REV.0104 - 0105 - 0106

Size: XS - S - M - L - XL - XXL

Color:

Composition: 74% polyamide,

26% elastan

Washing: Was












REV.0105



REV.0106













Designed by surgeons

Best seller

Size	Waist	Hips	Thigh	Mid Thigh
XS	47 - 59	76 - 92	40 - 50	34 - 42
S	57 - 69	84 - 100	47 - 57	38 - 46
М	67 - 79	92 - 108	54 - 64	42 - 50
L	77 - 89	100 - 116	61 - 71	46 - 54
XL	87 - 99	108 - 124	68 - 78	50 - 58
XXL	97 - 109	116 - 132	75 - 85	54 - 62
3XL	107 - 119	124 - 140	82 - 92	58 - 66

*NOTES: MEASUREMENTS ARE INDICATIVE AND EXPRESSED IN CENTIMETERS




LIPOSUCTION GARMENTS

Lipo Revée® is a collection of multifunctional garments that meet all the specific requirements of an effective compression treatment, after liposuction interventions in one or more areas of the body: abdominoplasty, thighplasty, torsoplasty, body reshaping, lipofilling, and buttock lift.

Thanks to the **Lycra® fabric**, these garments are able to self-mold ensuring the best fit and promoting the healing process.



Special liposuction support in the upper back







LIPO REVÉE

REV.0404



- 1 Triple front hook closure
- **2** Adjustable and removable straps
- 3 Wide hygienic opening
- 4 Abdominal support
- 5 Additional back support in Lycra®

Product code: REV.0404 - 0405 - 0406

Size: XS - S - M - L - XL - XXL - 3XL

Colors:

Composition: 74% polyamide,

26% elastan

Washing: $\[\] \boxtimes \[\] \boxtimes \[\] \boxtimes \[\] \otimes$





REV.0405



REV.0406













Best seller


Size	Waist	Hips	Thigh	Mid Thigh
XS	47 - 59	76 - 92	40 - 50	34 - 42
S	57 - 69	84 - 100	47 - 57	38 - 46
М	67 - 79	92 - 108	54 - 64	42 - 50
L	77 - 89	100 - 116	61 - 71	46 - 54
XL	87 - 99	108 - 124	68 - 78	50 - 58
XXL	97 - 109	116 - 132	75 - 85	54 - 62
3XL	107 - 119	124 - 140	82 - 92	58 - 66

*NOTES: MEASUREMENTS ARE INDICATIVE AND EXPRESSED IN CENTIMETERS



BBL | GARMENTS |

STANDARD HOLES



REV.0107 REV.0108








Lipo Revée® is a collection of multifunctional garments that meet all the specific requirements of an effective compression treatment, after liposuction interventions in one or more areas of the body: abdominoplasty, thighplasty, torsoplasty, body reshaping, lipofilling, and buttock lift. Thanks to the Lycra® fabric, these garments are able to self-mold ensuring the best fit and promoting the healing process.

BUTTOCKS&HIPS HOLES



REV.0407 REV.0408

REVÉE BBL



REVÉE BBL



1 Adjustable and removable straps

2 Abdominal support

3 Closure: 3 positions of hooks

4 Opening on the buttocks

5 Wide hygienic opening


Product code: REV.0107 - 0108 Size: XS - S - M - L - XL - XXL - 3XL

Colors:

Composition: 74% polyamide,

26% elastan

Washing: $\[\] \boxtimes \[\] \boxtimes \[\] \boxtimes \[\] \otimes$



STANDARD HOLES

REV.0107



REV.0108





compression

3XL



hooks



by surgeons

107 - 119







XS 47 - 59 76 - 92 40 - 50 34 - 42 S 57 - 69 84 - 100 47 - 57 38 - 46 M 67 - 79 92 - 108 54 - 64 42 - 50 L 77 - 89 100 - 116 46 - 54 61 - 71 68 - 78 50 - 58 XL 87 - 99 108 - 124 XXL 97 - 109 116 - 132 75 - 85 54 - 62

124 - 140

*NOTES: MEASUREMENTS ARE INDICATIVE AND EXPRESSED IN CENTIMETERS

82 - 92

58 - 66



REVÉE BBL



- 1 Adjustable and removable straps
- 2 Abdominal support
- 3 Triple front hook closure
- 4 Bigger buttonks
- 5 Wide hygienic opening
- 6 Additional back support in Lycra[®]

Product code: REV.0407 - 0408 Size: XS - S - M - L - XL - XXL - 3XL

Colors:

Composition: 74% polyamide,

26% elastan

Washing: $\[\[\] \] \[\] \[\$







BUTTOCKS&HIPS HOLES

REV.0407



REV.0408







hooks









Size	Waist	Hips	Thigh	Mid Thigh
XS	47 - 59	76 - 92	40 - 50	34 - 42
S	57 - 69	84 - 100	47 - 57	38 - 46
M	67 - 79	92 - 108	54 - 64	42 - 50
L	77 - 89	100 - 116	61 - 71	46 - 54
XL	87 - 99	108 - 124	68 - 78	50 - 58
XXL	97 - 109	116 - 132	75 - 85	54 - 62
3XL	107 - 119	124 - 140	82 - 92	58 - 66

TUMMY TUCK GARMENTS





After undergoing **abdominoplasty** surgery, it is essential to have an **support in the post-operative** period. In collaboration with our team of experts we have created a line of compression garments specially developed to improve your recovery after abdominoplasty surgery and emphasise the postoperative results.

Application fields: post-operative compression therapy following an operation of: hip liposuction abdominoplasty, torsoplasty and lipofilling with harvesting in the abdominal region, caesarean section.













BODY REVÉE

REV.1351



1 Adjustable and removable straps

2 Abdominal support

3 Side closure

Product code: REV.1351

Size: XXS - XS - S - M - L - XL - XXL - 3XL

Colors:

Composition: 74% polyamide,

26% elastan

Washing: ₩ ⋈ ⋈ ⋈ ⋈



Size	Waist	Hips
XXS	39 - 51	70 - 86
XS	49 - 61	78 - 94
S	59 - 71	86 - 102
М	69 - 81	94 - 110
L	79 - 91	102 - 118
XL	89 -101	110 - 126
XXL	99 -111	118 - 134
3XL	109 - 121	126 - 142











Designed Latex by surgeons free



Learn more







BODY REVÉE

REV.1352





2 Abdominal support

3 Front closure

4 Additional back support in Lycra®

Product code: REV.1352

Size: XS - S - M - L - XL - XXL - 3XL

Colors:

Composition: 74% polyamide,

26% elastan

Washing: 🖾 🔀 🖄 🛇



Size	Waist	Hips
XS	39 - 51	70 - 86
S	49 - 61	78 - 94
М	59 - 71	86 - 102
L	69 - 81	94 - 110
XL	79 - 91	102 - 118
XXL	89 -101	110 - 126
3XL	99 -111	118 - 134

REVEE

BODY REVÉE

REV.1353





2 Closure: 3 positions of hooks on both sides

3 Soft jacquard elastic bend by Revée®

Product code: REV.1353

Size: XS - S - M - L - XL - XXL - 3XL - 4XL

Colors:

Composition: 74% polyamide,

26% elastan

Washing: ₩ ⋈ ⋈ ⋈ ⋈

Use: post-surgical compression therapy after surgeries such as: caesarean section, hysterectomy, mini-abdominoplasty, lipofilling with extraction from the lower abdominal region, and lower-abdominal surgeries in general



Size	Hips
XS	74 - 80
S	81 - 87
М	88 - 94
L	95 - 101
XL	102 - 108
XXL	109 - 115
3XL	116 - 122
4XL	123 - 129







3 front hooks



Designed by surgeons post-parto





Regolabile sui Hips

















Diastasis recti



Adjustable velcro closure



Best seller







$R E V \check{E} E^{\circ}$

REVÉE BINDER

REV.0203 cm26





1 Breathable fabric

2 Contoured

3 Easy to put on

4 Abdominal support

5 Front cotton Rireinforcement

Product code: REV.0203 - 0204 Size: S - M - L - XL - XXL - 3XL

Composition: 87% polyamide,

13% elastan

Size	Hips
S	70 - 82
M	82 - 94
L	94 - 105
XL	105 - 117
XXL	117 - 130
3XL	130 - 150



















REVÉE SLEEVES







2 rows of hooks



Adjustable velcro closure



Best seller



Designed by surgeons



Learn more





REVÉE SLEEVES

REV.200





2 Closure: 2 positions of hooks

Product code: REV.200 Size: XS - S - M - L - XL - XXL

Colors:

Composition: 74% polyamide,

26% elastan

Washing: $\[\] \boxtimes \[\] \boxtimes \[\] \boxtimes \[\] \otimes$



Size	Arm
XS	24 - 25.5
S	25.5 - 27
M	27 - 28.5
L	28.5 - 30
XL	30 - 34
XXL	≥ 35

^{*}NOTES: MEASUREMENTS ARE INDICATIVE AND EXPRESSED IN CENTIMETERS

REVEE®

REVÉE SLEEVES

REV.201





- 1 Upper back compression
- 2 4 rows of hooks for upper adjustment
- 3 Triple front hooks for abdominal adjustment

Product code: REV.201 Size: XS - S - M - L - XL - XXL

Colors:

Composition: 74% polyamide,

26% elastan

Washing: \(\overline{40'} \overline{\overline{1}} \ov

Application fields: post-surgical compression therapy after brachioplasty, arm lifting surgery, armpit surgery and back liposuction.

Size	Arm	Under breast
XS	24 - 26	64 - 69
S	26 - 27	69 - 74
M	27 - 29	74 - 79
L	29 - 30	79 - 84
XL	30 - 33	84 - 91
XXL	33- 36	91 - 102



hooks



seller



Designed by surgeons















Unisex



FACE GARMENTS

Face compression garments provide targeted compression for chin liposuction or other facial surgical procedures. They are ideal for face and neck compression after neck lifts or chin liposuction. If applied correctly, these garments will increase blood circulation around the surgical area where

the surgery has been performed, as well as reduce swelling, eliminate any harmful bodily fluids, and speed up the healing process. In case of a sudden movement (e.g. sneeze or cough), these facial garments will minimize the discomfort



REVÉE MASKS

TOP OF HEAD-CHIN REV.5001



Product code: REV.5001

Size: S - M - L Colors:

Composition: 74% polyamide,

26% elastan

Washing: \bigcirc \bigcirc \bigcirc \bigcirc \bigcirc \bigcirc

Use: Lypoplasty of the chin, face lift.

Size	Top of head-chin
S	53 - 57
М	58 - 62
L	63 - 67







REV.5002



Product code: REV.5002

Size: S - M - L - XL

Colors:

Composition: 74% polyamide,

26% elastan

Washing: \bigcirc \bigcirc \bigcirc \bigcirc \bigcirc \bigcirc

Use: Forehead and face lift, lypoplasty of the chin, otoplasty, malar augmentation

Size	Top of head-chin
S	53 - 57
M	58 - 62
L	63 - 67
XL	68 - 72

REV.5003



Product code: REV.5003 Size: One-size-fits-all Colors:

Composition: 86% polyester,

14% rubber

Use: Lypoplasty of the chin, face lift, otoplasty.





Unisex



Abdominal compression Best seller



Cotton



Designed by surgeons

f O You

REVÉE PAD

Revée Pads are worn under the garments to provide **greater compression on the treated area** following reconstructive, aesthetic or reconstructive surgery. Product available in three different formats with different intended uses. The pads provide a surface smooth and soft surface in contact with the skin, improving comfort after liposuction or after a remodelling surgery body

contouring, help minimise the risk of haematomas, preventing the formation of bruises and skin folds, optimising lymphatic drainage by providing adequate compression accelerating and helping preventing healing water by retention and reducing the accumulation of fluid in stubborn areas. Also indicated for patients suffering from lipedema.





REVÉE PAD

REVÉE PAD01



Product code: REVÉE PAD01 Size: S - M / L - XL / XXL - 3XL

Colors:

Composition: 44% cotton,

56% polyurethane

Washing: ₩ 🔀 🖾 🖄 🛇

Size	GiroWaist
S - M	63 - 76
L - XL	77 - 94
XXL - 3XL	95 - 112





REVÉE PAD02



Product code: REVÉE PAD02

Size: One-size-fits-all

Colors:

Composition: 44% cotton,

56% polyurethane

Washing: $307 \bowtie 307 \bowtie 307 \bowtie 307$

REVÉE PAD03



Product code: REVÉE PAD03

Size: One-size-fits-all

Colors:

Composition: 44% cotton,

56% polyurethane

Washing: $\mathbf{X} \boxtimes \mathbf{X} \boxtimes \mathbf{X}$



2 pieces REVÉE PADO2 in each pack




THROMBOSAFE®













C	Circonferenza cm			Lunghezza cm	
J	сВ	сD	сG	iD	iG
Short	19 - 22	28 - 34	45 - 55	34 - 39	65 -75
Standard	19 - 22	28 - 34	45 - 55	39 - 44	75 - 85
Long	19 - 22	28 - 34	45 - 55	44 - 49	85 - 95

NA	Circonferenza cm			Lunghezza cm	
IVI	сВ	cD	сG	iD	iG
Short	22 - 25	34 - 40	50 - 60	34 - 39	65 -75
Standard	22 - 25	34 - 40	50 - 60	39 - 44	75 - 85
Long	22 - 25	34 - 40	50 - 60	44 - 49	85 - 95

		Circonferenza cm			Lunghezza cm	
L		сВ	сD	сG	iD	iG
Short		25 - 28	40 - 46	55 - 65	34 - 39	65 -75
Standard		25 - 28	40 - 46	55 - 65	39 - 44	75 - 85
Long		25 - 28	40 - 46	55 - 65	44 - 49	85 - 95

ΧI	Circonferenza cm			Lunghezza cm	
	сВ	cD	сG	iD	iG
Short	28 - 32	46 - 53	60 - 75	34 - 39	65 -75
Standard	28 - 32	46 - 53	60 - 75	39 - 44	75 - 85
Long	28 - 32	46 - 53	60 - 75	44 - 49	85 - 95

GYNECOMASTIA GARMENTS

Our Revée Male Collection provides a complete solution for post-operative recovery following gynecomastia, lipoplasty of the back, midriff and hips, liposuction and abdominal surgery. It combines an exclusive design, engineered with the support of world reknown plastic surgeons with the highest-

standards manufacturing process. Thanks to the unique characteristics of the Lycra[®] fabric (Latex free) it ensures optimal compression while providing maximum comfort to facillitate healing and repair of the tissues and inframammary structures.



REV.030



REV.032










REVÉE CHEST

REV.030





2 Pectoral support

3 Adjustable velcro closure


Product code: REV.030 Size: XXS - XS - S - M - L - XL

Colors:

Composition: 74% polyamide,

26% elastan

Washing: $\[\] \boxtimes \[\] \boxtimes \[\] \boxtimes \[\] \otimes$



Size	Torax
XXS	72 - 80
XS	80 - 88
S	89 - 97
M	98 - 106
L	107 - 115
XL	116 - 124

















Pectoral compression



Adjustable velcro closure



3 front hooks



Best seller

REVÉE CHEST

REV.031





2 Abdominal support

3 Pectoral support

4 Additional X-shape support Velcro closure

Product code: REV.031 Size: XXS - XS - S - M - L - XL

Colors:

Composition: 74% polyamide,

26% elastan

Washing: Washing: Washing: Washing



Size	Torax
XXS	72 - 80
XS	80 - 88
S	89 - 97
M	98 - 106
L	107 - 115
XL	116 - 124















REVÉE CARDIO

REV.032



Revée Cardio can aid recovery after open heart surgery or other chest surgeries. It reduces pain, decreases the risk of complications, minimizes friction between the two halves of the sternum, while allowing to breathe more easily.

- 1 4 adjustable velcro closures
- 2 Sternal support
- 3 Contoured binder
- 4 Adjustable and removable straps



Product code: REV.032 Size: S - M - L - XL - XXL

Colors:

Composition: 74% polyamide,

26% elastan

Size	Torax
S	78 - 88
M	89 - 98
L	98 - 108
XL	109 - 118
XXL	119 - 128







LIPO REVÉE FOR MAN

REV.0104





Size	Waist	Hips	Thigh	Mid Thigh
XS	47 - 59	76 - 92	40 - 50	34 - 42
S	57 - 69	84 - 100	47 - 57	38 - 46
М	67 - 79	92 - 108	54 - 64	42 - 50
L	77 - 89	100 - 116	61 - 71	46 - 54
XL	87 - 99	108 - 124	68 - 78	50 - 58
XXL	97 - 109	116 - 132	75 - 85	54 - 62
3XL	107 - 119	124 - 140	82 - 92	58 - 66

^{*}NOTES: MEASUREMENTS ARE INDICATIVE AND EXPRESSED IN CENTIMETERS

f S You





Our Revée Male Collection

provides a complete solution for post-operative recovery following gynecomastia, lipoplasty of the back, midriff and hips, liposuction and abdominal surgery. It combines an exclusive design, engineered with the support of world reknown plastic surgeons with the highest-standards manufacturing process. Thanks to the unique characteristics of the Lycra® fabric (Latex free) it ensures optimal compression while providing maximum comfort to facillitate healing and repair of the tissues and inframammary structures.

1 Closure: 3 positions of hooks

2 Adjustable and removable straps

3 Abdominal support

Product code: REV.0104 Size: XS - S - M - L - XL - XXL

Colors:

Composition: 74% polyamide,

26% elastan

Use: Abdominoplasty, liposuction of the abdomen, hips, thighs and lipofilling.



Designed

by surgeons Abdominal

compression

LIPO REVÉE FOR MAN

REV.0504





Size	Torax	Waist	Hips	Thigh
XS	87 - 91	75 - 79	82 - 85	47 - 50
S	92 - 96	80 - 86	86 - 90	51 - 55
М	97 - 102	87 - 93	91 - 97	56 - 60
L	103 - 109	94 - 101	98 - 104	61 - 65
XL	110 - 117	102 - 110	105 - 112	66- 70
XXL	118 - 125	111 - 119	113 - 120	71 - 75

^{*}NOTES: MEASUREMENTS ARE INDICATIVE AND EXPRESSED IN CENTIMETERS

f S You















Adjustable velcro closure



Best seller







REVÉE BINDER

REV.0203 0 cm26



- Breathable fabric
- 2 Contoured
- 3 Easy to put on
- 4 Abdominal support
- 5 Front cotton reinforcement

Product code: REV.0203 - 0204 Size: S - M - L - XL - XXL - 3XL

Colors:

Height: 26 cm - 32 cm

Composition: 87% polyamide,

13% elastan

Washing: 🕷 🛛 🖄





REV.0204 9cm32



Size	Hips
S	70 - 82
М	82 - 94
L	94 - 105
XL	105 - 117
XXL	117 - 130
3XL	130 - 150

*NOTES: MEASUREMENTS ARE INDICATIVE AND EXPRESSED IN CENTIMETERS







Medical silicone Odourless Without parabens



Silicone gels are the first treatment recommended by international guidelines to cure and prevent hypertrophic scars and keloids. (4)*



In the best pharmacies and healthcare shops





REVEE® SCARGEL

REVEE® SCAR GEL

TREATMENT FOR SCARS RECOMMENDED BY PLASTIC SURGEONS



It reduces redness and itching, while improving the aesthetical aspect. The scar becomes smoother and more elastic.



Odourless.



Quick-drying silicone protective film.

It is indicated for the treatment and prevention of hypertrophic or keloid scars caused by a burn, a surgical procedure or o trauma.


It is also used for immature and red scars in relief during the active healing phase.



How to use:

- 1 Clean the area to be treated. Then dry by dubbing;
- 2 Apply a thin layer of Revée Scar Gel. Do not massage;
- 3 After Revée[®] Scar Gel is dried you can apply: makeup, cosmetics or sunscreen;
- 4 Apply twice a day, morning and evening for at least 2/3 months.

Ingredients: Cyclopentasiloxane, Cyclohexasiloxane, dimethicone/vinyl dimethicone crosspolymer.









PROTECT SCARS FROM THE SUN.

GLUTEN FREE

www.revee.it





Sunscreen SPF 50+.

DERMATOLOGICALLY TESTED.

Highly resistant to water. High protection



Revée proudly supports research through the project

SUTURE IN SPACE

#SUTUREINSPACE



Video YouTube



Learn more







Azienda Ospedaliero Universitaria Careggi





UNIVERSITÀ DEGLI STUDI FIRENZE











Agenzia Spaziale Italiana





UNIVERSITÀ DEGLI STUDI FIRENZE

ACTA





















IN SPACE





Cute Project Onlus was founded in 2012 with the aim of offering theoretical and practical training to health personnel in developing countries, in the field of plastic reconstructive surgery, with specific address for the treatment of burns and their outcomes.

In recent years he has carried out several training missions in Africa, Congo, Benin, Uganda and started collaborative projects with the Serming of Turin and several Piedmont schools.

- It carries out humanitarian training missions for doctors and other local health figures through theoretical and practical lessons in the operating theater.
- Establishes scholarships in Italy, aimed at health personnel in developing countries and aimed at completing the process of acquiring specialist skills.
- It carries out health and educational activities for people in difficulty on the national and regional territory.

#REVEESUPPORT





Fiscal Code **97573900012**

PRIMARY PREVENTION

SECONDARY PREVENTION

TERTIARY PREVENTION

Interventions in schools, the fight against obesity, smoking, alcohol and the importance of self-examination. Risk factors for breast cancer.

Preventive examinations and guidance in-depth diagnostic.

Post-therapy social reintegration and psychological support.







Sede di Torino



#REVEESUPPORT









The Diastasis Women's OdV Association arises from the common effort of a group of mothers affected by diastasis of the abdominal rectum. The purpose of the Association is to inform as much as possible about this pathology, which is unfortunately still unknown today. Through events, collaborations, meetings, medical conferences, wants to give importance of women's health. Through the site and social accounts, they support women on their journey, from diagnosis / self-diagnosis to post-surgical, divulging practical and correct

information, which they make available thanks to doctors and professionals who work for free on their side. They organize awareness campaigns to recognize the pathology and treatment by the National Health Service in all Italian Regions. They provide assistance, psychological support and guidance to women who ask for help. They help women in post-partum, in choosing the most appropriate treatments for diastasis and for strengthening exercises, creating cooperation with professionals from all over Italy.





Information channel on plastic, reconstructive and aesthetic surgery









www.revee.news

REFERENCES

- 1 Chen et al. 1998 "The effect of a pressure garment on post-surgical drainage and seroma formation in breast cancer patients." Singapore Med J. 1998 Sep;39(9):412-5.
- 2 Rothman et al. 2014 "Abdominal binders may reduce pain and improve physical function after major abdominal surgery - a systematic review." Dan Med J. Nov;61(11):A4941.
- 3 Klement et al. 2010 "Supportive treatment using a compression garment vest of painful sternal instability following deep surgical wound infection: a case report." J Med Case Rep. 2010 Aug 11;4:266. doi: 10.1186/1752-1947-4-266.
- 4 Hsu et al. 2017 "Review of Silicone Gel Sheeting and Silicone Gel for the Prevention of Hypertrophic Scars and Keloids." Wounds. 2017 May; 29(5):154-158.



Note		





REVEE srl C.so Quintino Sella 131, 10132, Turin - Italy

DEALER